

insertion of a catheter or other medical instruments.

The insertion of urinary catheters by the helps of Dispogel Urethral Lubricant Gel with Lidocaine HCL (12.5g and 8.5g) health professionals minimise associated trauma, pain and catheter-associated infections. After instillation of gel - WAIT 5-10 minutes prior to insertion of a catheter or other medical instruments time needed for achieving the lidocaine anesthetic effect.

Dispogel is a Class III medical device.

8.5g	12.5g	STEF	RILE	R (€ 2292
GMDN Code	Product Description		Brand	
37717	12.5g Accordion Bottle		Dis	pogel 12.5g
37717	8.5g Accordion Bottle		Di	spogel 8.5g







Dispogel 8.5g (For Female)



## DİSPOFARMA İLAÇ SANAYİ VE TİCARET LTD. ŞTİ.

İvedik Organize Sanayi Belgesi, Ağaç İşleri Sitesi, 30. Cad. No: 51 Yenimahalle - ANKARA / TURKEY Tel: +90 312 394 45 52 Fax: +90 312 394 45 58 info@dispofarma.com - www.dispofarma.com