

Method of Use
For Male



Step 1



Step 2



Step 3



Step 4



Step 5



Step 6



Step 7



Step 8

For Female



Step 1



Step 2



Step 3



Step 4



Step 5

Dispogel®



Dispogel®
Urethral Lubricant Gel with
2% Lidocaine HCL
11 ml & 6 ml

Dispogel® is a sterile lubricant gel containing local anesthetic 2% Lidocaine HCL and antiseptic (0.05 % Clorhexidine Digluconate). Dispogel is available in both 6ml and 11 ml single use pre-filled syringes each presented in a sterile package.

Dispogel® is used for instillation into the urethra prior to insertion of a catheter or other medical instruments. The insertion of urinary catheters by the helps of Dispogel® Urethral Lubricant Gel with 2% Lidocaine HCL (11 ml & 6 ml) health professionals minimise associated trauma, pain and catheter-associated infections.

After instillation of gel - WAIT 5-10 minutes prior to insertion of a catheter or other medical instruments- time needed for achieving the lidocaine anesthetic effect.

Dispogel® is a Class III medical device.



6ml	11ml	STERILE R	CE 2292
GMDN Code	Product Description	Brand	Pack Size
37717	6 ml Syringe	Dispogel®	Pack of 50
37717	11 ml Syringe	Dispogel®	Pack of 25



Dispofarma İlaç Sanayi ve Ticaret Ltd. Şti.
İvedik Organize Sanayi Bölgesi, Ağaç İşleri Sitesi, 1333. Cad. No: 51
Yenimahalle - ANKARA / TURKEY
Tel: +90 (0) 312 394 45 52 Fax:+90 (0) 312 394 45 58
info@dispofarma.com
www.dispofarma.com.tr

ANAESTHETIC

with

ANTISEPTIC

with

STERILE LUBRICANT