



Dispogel[®]

**Urethral Lubricant Gel
with 2 % Lidocaine HCL**

Dispogel[®] is a sterile lubricat gel containing local anesthetic (2% Lidocaine HCL and antiseptic (0.05 % Clorhexidine Digluconate). Dispogel[®] is available in 12.5 gr and 8,5 gr single use accordion bottle presented in a sterile package.

Dispogel[®] is ustü for instillation into the urethra prior to insertion of a catheter or other medical instruments. The insertion of urinary catheters by the helps of Dispogel[®]

Urethral Lubricant Gel with 2% Lidocaine HCL (12.5 gr and 8,5 gr) health professionals minimise associated trauma, pain and catheter-associated infections.

After instillation of gel - WAIT 5-10 minutes prior to insertion of a catheter or other medical instruments- time needed for achieving the lidocaine anesthetic effect.

Dispogel[®] is a Class III medical device.



12.5 gr



8.5 gr

8.5 gr	12.5 gr	STERIL R	CE 2292
GMDN Code	Product Description	Brand	Pack Size
37717	8.5 gr Accordion Bottle	Dispogel [®]	Pack of 100
37717	12.5 gr Accordion Bottle	Dispogel [®]	Pack of 100

Method of Use
For Male



Step 1



Step 2



Step 3



Step 4




Step 5



Step 6



Step 7



Step 8

ANAESTHETIC

with

ANTISEPTIC

with

STERILE LUBRICANT



Dispofarma İlaç Sanayi ve Ticaret Ltd. Şti.

İvedik Organize Sanayi Bölgesi, Ağaç İşleri Sitesi, 1333. Cad. No: 51
Yenimahalle - ANKARA / TURKEY

Tel: +90 (0) 312 394 45 52 Fax: +90 (0) 312 394 45 58

info@dispofarma.com

www.dispofarma.com.tr